



**SHOWN ZBU1PK**

Part#	Description	Bore		L		L2		L3		L4		D		X1		X2		X3
		mm	Inch	mm	Inch	mm	Inch	mm	Inch	mm	Inch	mm	Inch	mm	Inch	mm	Inch	
ZBU.5FPK	Union 1/32", fingertight	Ø0.25	Ø.010"	-32.5	-1.28"	17.53	.690"	11.18	.440"	4.77	.188"	Ø8	Ø.3125"	Ø11.1	Ø.437"	Ø6.35	Ø.25"	1/2" Hex
ZBU.5LFPK	Union 1/32", fingertight	Ø0.50	Ø.020"	-32.5	-1.28"	17.53	.690"	11.18	.440"	4.77	.188"	Ø8	Ø.3125"	Ø11.1	Ø.437"	Ø6.35	Ø.25"	1/2" Hex
ZBU.5TFPK	Union 1/32", fingertight	Ø0.79	Ø1/32"	-31.8	-1.25"	16.76	.660"	10.41	.410"	4.77	.188"	Ø8	Ø.3125"	Ø11.1	Ø.437"	Ø6.35	Ø.25"	1/2" Hex
ZBU1CFPK	Union 1/16", fingertight	Ø0.25	Ø.010"	-60.8	-2.39"	22.1	.870"	12.57	.495"	3.18	.125"	Ø9.6	Ø.3750"	11.1	.437"	Ø9.65	Ø.38"	1/2" Hex
ZBU1MFPK	Union 1/16", fingertight	Ø0.50	Ø.020"	-60.8	-2.39"	22.1	.870"	12.57	.495"	3.18	.125"	Ø9.6	Ø.3750"	11.1	.437"	Ø9.65	Ø.38"	1/2" Hex
ZBU1FPK	Union 1/16", fingertight	Ø0.75	Ø.030"	-60.8	-2.39"	22.1	.870"	12.57	.495"	3.18	.125"	Ø9.6	Ø.3750"	11.1	.437"	Ø9.65	Ø.38"	1/2" Hex
ZBU1TFPK	Union 1/16", fingertight	Ø1.59	Ø1/16"	-59.4	-2.34"	20.83	.820"	11.30	.445"	3.18	.125"	Ø9.6	Ø.3750"	11.1	.437"	Ø9.65	Ø.38"	1/2" Hex
ZBU1CPK	Union 1/16", hex	Ø0.25	Ø.010"	-38.0	-1.50"	22.1	.870"	12.57	.495"	3.18	.125"	Ø9.6	Ø.3750"	11.1	.437"	6.35 Hex	1/4" Hex	1/2" Hex
ZBU1MPK	Union 1/16", hex	Ø0.50	Ø.020"	-38.0	-1.50"	22.1	.870"	12.57	.495"	3.18	.125"	Ø9.6	Ø.3750"	11.1	.437"	6.35 Hex	1/4" Hex	1/2" Hex
ZBU1PK	Union 1/16", hex	Ø0.75	Ø.030"	-38.0	-1.50"	22.1	.870"	12.57	.495"	3.18	.125"	Ø9.6	Ø.3750"	11.1	.437"	6.35 Hex	1/4" Hex	1/2" Hex
ZBU1TPK	Union 1/16", hex	Ø1.59	Ø1/16"	-36.8	-1.45"	20.83	.820"	11.30	.445"	3.18	.125"	Ø9.6	Ø.3750"	11.1	.437"	6.35 Hex	1/4" Hex	1/2" Hex
ZBU2PK	Union 1/8", hex	Ø0.75	Ø.030"	-53.2	-2.10"	33.66	1.325"	22.86	.900"	6.35	.250"	Ø11.1	Ø.4375"	12.7	.500"	9.53 Hex	3/8" Hex	11/16" Hex
ZBU2MPK	Union 1/8", hex	Ø1.50	Ø.060"	-53.2	-2.10"	33.66	1.325"	22.86	.900"	6.35	.250"	Ø11.1	Ø.4375"	12.7	.500"	9.53 Hex	3/8" Hex	11/16" Hex
ZBU2LPK	Union 1/8", hex	Ø2.00	Ø.080"	-53.2	-2.10"	33.66	1.325"	22.86	.900"	6.35	.250"	Ø11.1	Ø.4375"	12.7	.500"	9.53 Hex	3/8" Hex	11/16" Hex
ZBU2TPK	Union 1/8", hex	Ø3.18	Ø1/8"	-51.2	-2.02"	31.75	1.250"	20.96	.825"	6.35	.250"	Ø11.1	Ø.4375"	12.7	.500"	9.53 Hex	3/8" Hex	11/16" Hex

THIS DOCUMENT AND THE INFORMATION WHICH IT CONTAINS SHALL NOT BE USED, EXPLOITED OR SOLD, AND SHALL NOT BE REPRODUCED OR DISCLOSED TO OTHERS WITHOUT THE EXPRESS WRITTEN PERMISSION OF VICI. THIS DOCUMENT SHALL REMAIN THE PROPERTY OF VICI AND SHALL BE RETURNED UPON DEMAND.

Änderungen / Modifications: a) 21.01.09 FP #568, added information

Material: KUNDE

Production: EN 602268-1 ISO 2768

Scale: 5:1

Customer Order: 29.02.2008 F.P.

Drawing No.: 21.01.2009 F.C.o

File No.: 506326.idw

Customer: KUNDE

Drawing No.: 1 - GE - 5063.26 a

**Bulkhead internal unions high pressure PEEK**

General dimensions

VICI VICI AG International, CH-6214 Schenkon