



A-A



Shown ZNF.5FPK

B-B



Shown ZNF1PKG

Dimension in mm [Inch]

THIS DOCUMENT AND THE INFORMATION WHICH IT CONTAINS SHALL NOT BE USED, EXPLOITED OR SOLD, AND SHALL NOT BE REVEALED OR DISCLOSED TO OTHERS WITHOUT THE EXPRESS WRITTEN PERMISSION OF VICI. THIS DOCUMENT SHALL REMAIN THE PROPERTY OF VICI AND SHALL BE RETURNED UPON DEMAND.

Description	Glass-filled PEEK ferrule Prod No	PEEK ferrule Prod No	CTFE ferrule Prod No	L	X
1/32" fingertight	ZNF.5FPKG	ZNF.5FPK		14.5 [.57"]	Ø6.35 [Ø.25"]
1/16" fingertight	ZNF1FPKG	ZNF1FPK	ZNF1FKF	26.9 [1.06"]	Ø9.4 [Ø.37"]
1/16" Hex, short	ZNF1PKG	ZNF1PK	ZNF1KF	16.3 [.64"]	6.35 [1/4"]
1/16" Hex, medium	MZNF1PKG	MZNF1PK	MZNF1KF	20.8 [.82"]	6.35 [1/4"]
1/16" Hex, long	LZNF1PKG	LZNF1PK	LZNF1KF	27.2 [1.07"]	6.35 [1/4"]

Aenderungen:
MODIFICATIONS:
PART NO.:

Material MATERIAL	Dimensionen DIMENSIONS
Projektion PROJECTION	Allgemein Toleranzen GENERAL TOLERANCES
SN EN 22768-1/ ISO 2768	

No twist one-piece fittings for 1/32" and 1/16" tubing
General dimensions

Masstab SCALE	Gezeichnet DRAWN	03.10.2012	F.P.
4:1	Geprüft / Freigabe CHECKED/ RELEASE	03.10.2012	FCo
Unterverz.: SUB DIREC.:		KUNDE	
Zeichnungs Nr. / DRAWING NO.		3 - GE - 5063.96	

VICI[®] VICI AG International, CH-6214 Schenkon