


C42-9786A valve, 6 port 2-pos., shown

Dimension in mm [Inch]



1/2-20 UNF threaded fitting details for use with 1/4" OD tubing.

Part#	Description	Bore		X		Y		L	
		mm	Inch	mm	Inch	mm	Inch	mm	Inch
C42-9784A	Valve-plug, 4 port	Ø4.57	Ø0.18"	Ø37.85	Ø1.49"	-70	-2.76"	48	1.89"
C42-9786A	Valve-plug, 6 port	Ø4.57	Ø0.18"	Ø50.55	Ø1.99"	-83	-3.26"	51.3	2.02"

THIS DOCUMENT AND THE INFORMATION WHICH IT CONTAINS SHALL NOT BE USED, EXPLOITED OR SOLD, AND SHALL NOT BE REVEALED OR DISCLOSED TO OTHERS WITHOUT THE EXPRESS WRITTEN PERMISSION OF VICI. THIS DOCUMENT SHALL REMAIN THE PROPERTY OF VICI AND SHALL BE RETURNED UPON DEMAND.

Aenderungen:
MODIFICATIONS:

PART NO.:

Material
MATERIAL

Dimensionen
DIMENSIONS

Projektion
PROJECTION

Allgemein Toleranzen
GENERAL TOLERANCES

SN EN 22768-1/ISO 2768

**C42 Valve, 1/2-20 UNF
Air actuated**

General dimensions

Masstab SCALE	Gezeichnet DRAWN	23.10.2009	F.P.
1:1	Gepflicht / Freigabe CHECKED/RELEASE	23.10.2009	FCo
CAD Nr. CAD NO.	506353.idw		

Unterverz.:
SUB DIRC.: KUNDE

Zeichnungs Nr. / DRAWING NO.



VICI AG International, CH-6214 Schenkon

3 - GE - 5063.53