



Part Nr.	Description	X	L
C25G-24520EMT	20 Position	Ø44.45 mm [Ø1.75"]	~41.7 mm [~1.64"]
C25G-24524EMT	24 Position	Ø44.45 mm [Ø1.75"]	~41.7 mm [~1.64"]
C25G-24528EMT	28 Position	Ø44.45 mm [Ø1.75"]	~41.7 mm [~1.64"]

THIS DOCUMENT AND THE INFORMATION WHICH IT CONTAINS SHALL NOT BE USED, EXPLOITED OR SOLD, AND SHALL NOT BE REVEALED OR DISCLOSED TO OTHERS WITHOUT THE EXPRESS WRITTEN PERMISSION OF VICI. THIS DOCUMENT SHALL REMAIN THE PROPERTY OF VICI AND SHALL BE RETURNED UPON DEMAND.

Aenderungen: MODIFICATIONS:

PART NO.: C25G-2452*EMT

Model C25G, low pressure Stream selector, 1/16" fittings

with micro electric actuator. General dimensions

VICI® VICI AG International, CH-6214 Schenkon

Masstab SCALE	Gezeichnet DRAWN	08.05.2013	F.P.
1:1	Geprüft / Freigabe CHECKED/ RELEASE	08.05.2013	FCo
	CAD Nr. CAD NO.	506406.idw	

Unterverz.: SUB DIREC.: KUNDE

Zeichnungs Nr. / DRAWING NO. 3 - GE - 5064.06

Material MATERIAL Dimensionen DIMENSIONS

Projektion PROJECTION Allgemein Toleranzen GENERAL TOLERANCES

SN EN 22768-1/ ISO 2768