




Dimension in mm [Inch]



Micro lynx controller Part# CP-CM1-P



Power supply: 24 VDC, 110/220 VAC, CE
Part# PS24VDC-CE
35.6 x 59.7 x 118 [1.4" x 2.35" x 4.65"]

THIS DOCUMENT AND THE INFORMATION WHICH IT CONTAINS SHALL NOT BE USED, EXPLOITED OR SOLD, AND SHALL NOT BE REVEALED OR DISCLOSED TO OTHERS WITHOUT THE EXPRESS WRITTEN PERMISSION OF VICI. THIS DOCUMENT SHALL REMAIN THE PROPERTY OF VICI AND SHALL BE RETURNED UPON DEMAND.		Material	Dimensionen
Anderungen: MODIFICATIONS:		PROJECTION	Allgemein Toleranzen
PART NO.: CP2-4841-100M1		CAD Nr.	GENERAL TOLERANCES
M6 Pump with controller and stepper motor		CAD NO.	SN EN 22768-1/ ISO 2768
General dimensions		Masstab SCALE	Gezeichnet DRAWN
VICI® VICI AG International, CH-6214 Schenkon		1:1	16.06.2009 F.P.
Unterverz.: SUB DIREC.:		KUNDE	Geprüft / Freigabe CHECKED/ RELEASE
Zeichnungs Nr. / DRAWING NO.		16.06.2009 FCo	
2 - GE - 5063.51		506351.idw	