



SHOWN EZRU21

Part#	Description	Bore		L		X1	X2	X3
		mm	Inch	mm	Inch			
EZRU1.5	External / internal reducing union 1/16" to 1/32"	Ø0.25	Ø.010"	-23.6	-.93"	1/4" Hex	1/4" Hex	3/16" Hex
EZRU1.5L	External / internal reducing union 1/16" to 1/32"	Ø0.50	Ø.020"	-23.6	-.93"	1/4" Hex	1/4" Hex	3/16" Hex
EZRU1.5T	External / internal reducing union 1/16" to 1/32"	Ø0.89	Ø.035"	-22.6	-.89"	1/4" Hex	1/4" Hex	3/16" Hex
EZRU2.5	External / internal reducing union 1/8" to 1/32"	Ø0.25	Ø.010"	-27.0	-1.06"	3/8" Hex	7/16" Hex	3/16" Hex
EZRU2.5L	External / internal reducing union 1/8" to 1/32"	Ø0.50	Ø.020"	-27.0	-1.06"	3/8" Hex	7/16" Hex	3/16" Hex
EZRU2.5T	External / internal reducing union 1/8" to 1/32"	Ø0.89	Ø.035"	-26.0	-1.03"	3/8" Hex	7/16" Hex	3/16" Hex
EZRU21C	External / internal reducing union 1/8" to 1/16"	Ø0.25	Ø.010"	-32.8	-1.29"	3/8" Hex	7/16" Hex	1/4" Hex
EZRU21	External / internal reducing union 1/8" to 1/16"	Ø0.75	Ø.030"	-33.3	-1.31"	3/8" Hex	7/16" Hex	1/4" Hex
EZRU21T	External / internal reducing union 1/8" to 1/16"	Ø1.65	Ø.065"	-32.0	-1.26"	3/8" Hex	7/16" Hex	1/4" Hex
EZRU4.5	External / internal reducing union 1/4" to 1/32"	Ø0.25	Ø.010"	-30.8	-1.21"	9/16" Hex	9/16" Hex	3/16" Hex
EZRU4.5T	External / internal reducing union 1/4" to 1/32"	Ø0.89	Ø.035"	-29.8	-1.17"	9/16" Hex	9/16" Hex	3/16" Hex
EZRU41C	External / internal reducing union 1/4" to 1/16"	Ø0.25	Ø.010"	-36.0	-1.42"	9/16" Hex	9/16" Hex	1/4" Hex
EZRU41	External / internal reducing union 1/4" to 1/16"	Ø0.75	Ø.030"	-36.5	-1.44"	9/16" Hex	9/16" Hex	1/4" Hex
EZRU41T	External / internal reducing union 1/4" to 1/16"	Ø1.65	Ø.065"	-35.3	-1.39"	9/16" Hex	9/16" Hex	1/4" Hex
EZRU42	External / internal reducing union 1/4" to 1/8"	Ø1.0	Ø.040"	-43.0	-1.69"	9/16" Hex	9/16" Hex	3/8" Hex
EZRU42L	External / internal reducing union 1/4" to 1/8"	Ø2.0	Ø.080"	-43.0	-1.69"	9/16" Hex	9/16" Hex	3/8" Hex
EZRU42T	External / internal reducing union 1/4" to 1/8"	Ø3.18	Ø.125"	-41.0	-1.62"	9/16" Hex	9/16" Hex	3/8" Hex

THIS DOCUMENT AND THE INFORMATION WHICH IT CONTAINS SHALL NOT BE USED, EXPLOITED OR SOLD, AND SHALL NOT BE REPRODUCED OR DISCLOSED TO OTHERS WITHOUT THE EXPRESS WRITTEN PERMISSION OF VICI. THIS DOCUMENT SHALL REMAIN THE PROPERTY OF VICI AND SHALL BE RETURNED UPON DEMAND.

Material	Dimensions
Position	General Tolerances
Part No.	EN 62278-1 ISO 2768
Material	Material
Scale	14.12.2005 F.P.
5:1	15.12.2005 F.C.
CAD No.	506055.idw
Customer	KUNDE
Drawn by	Zürcherstr. 100 / DRAWING NO.
VICI® VICI AG International, CH-6214 Schenkon	
1 - GE - 5060.55	