


**SHOWN ZRU21**



Part#	Description	Bore		L1		L2		L3		L4		X1	X2	X3
		mm	Inch	mm	Inch	mm	Inch	mm	Inch	mm	Inch			
ZRU1.5XC	1/16" to 1/32" reducing union	Ø0.15	Ø.006"	-31.4	-1.235"	19.18	0.755"	-8.0	-0.313"	-4.2	-0.167"	1/4" Hex	1/4" Hex	3/16" Hex
ZRU1.5	1/16" to 1/32" reducing union	Ø0.25	Ø.010"	-31.8	-1.250"	19.56	0.770"	-8.0	-0.313"	-4.2	-0.167"	1/4" Hex	1/4" Hex	3/16" Hex
ZRU1.5L	1/16" to 1/32" reducing union	Ø0.50	Ø.020"	-32.3	-1.270"	20.07	0.790"	-8.0	-0.313"	-4.2	-0.167"	1/4" Hex	1/4" Hex	3/16" Hex
ZRU1.5T	1/16" to 1/32" reducing union	Ø0.89	Ø.035"	-31.0	-1.220"	18.80	0.740"	-8.0	-0.313"	-4.2	-0.167"	1/4" Hex	1/4" Hex	3/16" Hex
ZRU2.5	1/8" to 1/32" reducing union	Ø0.25	Ø.010"	-37.7	-1.485"	23.75	0.935"	-9.7	-0.383"	-4.2	-0.167"	7/16" Hex	3/8" Hex	3/16" Hex
ZRU2.5L	1/8" to 1/32" reducing union	Ø0.50	Ø.020"	-38.2	-1.505"	24.26	0.955"	-9.7	-0.383"	-4.2	-0.167"	7/16" Hex	3/8" Hex	3/16" Hex
ZRU2.5T	1/8" to 1/32" reducing union	Ø0.89	Ø.035"	-37.0	-1.455"	23.0	0.905"	-9.7	-0.383"	-4.2	-0.167"	7/16" Hex	3/8" Hex	3/16" Hex
ZRU21C	1/8" to 1/16" reducing union	Ø0.25	Ø.010"	-43.5	-1.711"	25.78	1.015"	-9.7	-0.383"	-8.0	-0.313"	7/16" Hex	3/8" Hex	1/4" Hex
ZRU21	1/8" to 1/16" reducing union	Ø0.75	Ø.030"	-44.6	-1.756"	26.92	1.060"	-9.7	-0.383"	-8.0	-0.313"	7/16" Hex	3/8" Hex	1/4" Hex
ZRU21L	1/8" to 1/16" reducing union	Ø1.0	Ø.040"	-44.6	-1.756"	26.92	1.060"	-9.7	-0.383"	-8.0	-0.313"	7/16" Hex	3/8" Hex	1/4" Hex
ZRU21T	1/8" to 1/16" reducing union	Ø1.65	Ø.065"	-42.7	-1.681"	25.02	0.985"	-9.7	-0.383"	-8.0	-0.313"	7/16" Hex	3/8" Hex	1/4" Hex
ZRU41C	1/4" to 1/16" reducing union	Ø0.25	Ø.010"	-48.4	-1.903"	29.59	1.165"	-10.8	-0.425"	-8.0	-0.313"	9/16" Hex	1/2" Hex	1/4" Hex
ZRU41	1/4" to 1/16" reducing union	Ø0.75	Ø.030"	-49.5	-1.948"	30.73	1.210"	-10.8	-0.425"	-8.0	-0.313"	9/16" Hex	1/2" Hex	1/4" Hex
ZRU41T	1/4" to 1/16" reducing union	Ø1.65	Ø.065"	-47.6	-1.873"	28.83	1.135"	-10.8	-0.425"	-8.0	-0.313"	9/16" Hex	1/2" Hex	1/4" Hex
ZRU42	1/4" to 1/8" reducing union	Ø0.75	Ø.030"	-55.4	-2.183"	34.93	1.375"	-10.8	-0.425"	-9.7	-0.383"	9/16" Hex	1/2" Hex	3/8" Hex
ZRU42L	1/4" to 1/8" reducing union	Ø2.0	Ø.080"	-55.4	-2.183"	34.93	1.375"	-10.8	-0.425"	-9.7	-0.383"	9/16" Hex	1/2" Hex	3/8" Hex
ZRU42T	1/4" to 1/8" reducing union	Ø3.25	Ø.128"	-53.5	-2.108"	33.02	1.300"	-10.8	-0.425"	-9.7	-0.383"	9/16" Hex	1/2" Hex	3/8" Hex

THIS DOCUMENT AND THE INFORMATION WHICH IT CONTAINS SHALL NOT BE USED, EXPLOITED OR SOLD, AND SHALL NOT BE REPRODUCED OR DISCLOSED TO OTHERS WITHOUT THE EXPRESS WRITTEN PERMISSION OF VICI. THIS DOCUMENT SHALL REMAIN THE PROPERTY OF VICI AND SHALL BE RETURNED UPON DEMAND.

Material	Dimensions
Production	General Tolerances
20.03.2006 F.P. 27.03.2006 F.C.o. 506077 idw	EN EN 22768-1 ISO 2768 KUNDE

**Internal reducing unions**  
**Stainless steel**  
 General dimensions

**VICI** VICI AG International, CH-6214 Schenkon  
**1 - GE - 5060.77**