



SHOWN ZUFR1

Standard porosity 2 µm

Part#	Description	Bore		L		X1	X2	X3	X4	Replacement screen / frit
		mm	Inch	mm	Inch					
ZUFR.5	Filter union 1/32" with screen	∅0.25	∅.010"	~36.5	~1.440"	1/4" Hex	3/16" Hex	1/4" Hex	3/16" Hex	2SR1
ZUFR.5F	Filter union 1/32" with frit	∅0.25	∅.010"	~37.5	~1.475"	1/4" Hex	3/16" Hex	1/4" Hex	3/16" Hex	2FR1
ZUFR1	Filter union 1/16" with screen	∅0.50	∅.020"	~47.4	~1.865"	7/16" Hex	1/4" Hex	3/8" Hex	1/4" Hex	2SR2
ZUFR1F	Filter union 1/16" with frit	∅0.50	∅.020"	~48.7	~1.915"	7/16" Hex	1/4" Hex	3/8" Hex	1/4" Hex	2FR2
ZUFR1C	Filter union 1/16" with screen	∅0.25	∅.010"	~47.0	~1.850"	7/16" Hex	1/4" Hex	3/8" Hex	1/4" Hex	2SR2
ZUFR1CF	Filter union 1/16" with frit	∅0.25	∅.010"	~48.0	~1.890"	7/16" Hex	1/4" Hex	3/8" Hex	1/4" Hex	2FR2
ZRUFR1.5	Reducing filter union 1/16" to 1/32" with screen	∅0.25	∅.010"	~41.8	~1.645"	1/4" Hex	1/4" Hex	1/4" Hex	3/16" Hex	2SR1
ZRUFR1.5F	Reducing filter union 1/16" to 1/32" with frit	∅0.25	∅.010"	~42.7	~1.680"	1/4" Hex	1/4" Hex	1/4" Hex	3/16" Hex	2FR1
ZUFR2	Filter union 1/8" with screen	∅2.00	∅.080"	~67.3	~2.650"	9/16" Hex	3/8" Hex	9/16" Hex	3/8" Hex	2SR4
ZUFR2F	Filter union 1/8" with frit	∅2.00	∅.080"	~67.8	~2.670"	9/16" Hex	3/8" Hex	9/16" Hex	3/8" Hex	2FR4
ZRUFR21	Reducing filter union 1/8" to 1/16" with screen	∅0.75	∅.029"	~53.7	~2.115"	7/16" Hex	3/8" Hex	3/8" Hex	1/4" Hex	2SR2
ZRUFR21F	Reducing filter union 1/8" to 1/16" with frit	∅0.75	∅.029"	~54.7	~2.155"	7/16" Hex	3/8" Hex	3/8" Hex	1/4" Hex	2FR2
ZRUFR41	Reducing filter union 1/4" to 1/16" with screen	∅1.00	∅.039"	~61.4	~2.415"	9/16" Hex	1/2" Hex	1/2" Hex	1/4" Hex	2SR4
ZRUFR41F	Reducing filter union 1/4" to 1/16" with frit	∅1.00	∅.039"	~62.4	~2.455"	9/16" Hex	1/2" Hex	1/2" Hex	1/4" Hex	2FR4
ZRUFR42	Reducing filter union 1/4" to 1/8" with screen	∅2.00	∅.080"	~72.7	~2.860"	9/16" Hex	1/2" Hex	9/16" Hex	3/8" Hex	2SR4
ZRUFR42F	Reducing filter union 1/4" to 1/8" with frit	∅2.00	∅.080"	~73.7	~2.900"	9/16" Hex	1/2" Hex	9/16" Hex	3/8" Hex	2FR4

THIS DOCUMENT AND THE INFORMATION WHICH IT CONTAINS SHALL NOT BE USED, EXPLOITED OR SOLD, AND SHALL NOT BE REPRODUCED OR DISCLOSED TO OTHERS WITHOUT THE EXPRESS WRITTEN PERMISSION OF VICI. THIS DOCUMENT SHALL REMAIN THE PROPERTY OF VICI AND SHALL BE RETURNED UPON DEMAND.

Material:	Dimensions:
Production:	General Tolerances:
EN ISO 2768-V ISO 2768	

Part No.:

Scale: 5:1

Revision: 07.09.2006 F.C.O.

Drawn: 506090.idw

Customer: KUNDE

Drawing No. / Drawing No.:

VICI VICI AG International, CH-6214 Schenkon

1 - GE - 5060.90